

InnoVet Select



Specifications

- ~ kVp range: 40 - 125, 64 selections
- ~ mA stations: 100S, 300 L
- ~ Timer range: 1/120th - 1 sec., 18 selections
- ~ mAs range: 1.0 - 300 mAs
- ~ Digital mAs readout
- ~ Digital kVp readout
- ~ Digital diagnostic readout
- ~ Tilt-out control panel
- ~ Built-in tube protection
- ~ Solid state adjustable, audible exposure signal
- ~ SCR contactor
- ~ UL classified

Table/Tubestand

- ~ Enclosed, integrated radiographic table
- ~ Table top: 53" x 24"
- ~ Welded base construction
- ~ Transformer mounting shelf
- ~ Grid cabinet with Simplex™ cassette tray
- ~ Grid: 8:1 ratio, 103 lpi
- ~ Grid cabinet/tubestand interlock
- ~ Angulation dial and operator handgrips
- ~ Foot switch with 20 ft. cord
- ~ UL classified



X-Ray Tube

- ~ 1.0 - 2.0 mm F.S.
- ~ 140,000 H.U. with 90 degree horns

Cables

- ~ 10 ft. high voltage cables with federal terminations

Collimator

- ~ Manual collimator with light field
- ~ Swivel mount
- ~ Field light "on" time adjustable to 30 - 120 seconds
- ~ Long-life "soft-start" bulb circuits
- ~ Meets federal standards for radiation leakage protection
- ~ Easy to operate with gloves on