

VetTeKLC



Universal

VetTekLC



VetTekLC is the most powerful, rugged, feature enhanced and affordable veterinary system on the market. **VetTekLC** was designed to meet the needs of the practitioner on a limited budget without having to compromise features that ease the process of obtaining quality radiographs.

A robust table with floating tubestand provides benefits of accurate and effortless positioning of uncooperative patients. The added feature of our 400 mA table based generator allows use of shorter exposure times and setting techniques without leaving the patient.

Standard Features

- Stationary Top Table, 24" x 56"
- Table Based 400 mA Generator
- Digital mAs and kVp Displays on Tubestand
- 60-Second Light on Collimator
- Collimator Swivel, aids in performing split film table top views
- Collimator Tape Measure, accurately measures for table-top lateral or angled views
- Two Stage Footswitch
- Grid Cabinet with 103 line, 8:1 ratio grid and cassette tray. Accepts all size cassettes, both inch and metric
- Grid Cabinet/Tubestand Interlock, centers beam to cassette
- Adjustable Tube Height, provides for table top, in cabinet and table-top lateral views
- Detent for 40" Table Top and in Table SID (source, image, distance)
- Manual Locks for Tube and Tubestand
- Automatic Compensation of mA and kVp, no need to manually adjust for line voltage fluctuations
- Transformer Shelf, elevates transformer off floor

Technical Specifications

- Maximum mA: 400
- mA Stations: 4, 50S/100S and 300L/400L
- kVp Range: 40 to 125
- Time Stations: 13, 1/120th to 1 second

Made in U.S.A. – An ISO 13485 Certified Facility

Universal

Serving the Veterinary Profession for Over 75 years

www.delmedical.com



Medical
Imaging
Corporation

11550 W. King St. • Franklin Park, IL 60131 • Phone: 800-800-6006 • Fax: 847-288-7011